

Abstract

Background: Ceftaroline (CPT), the active form of CPT fosamil, has *in vitro* activity against resistant Gram-positive organisms, including MRSA and multidrug-resistant (R) *S. pneumoniae* (SPN). CPT was FDA-approved in late 2010 for the indications of community-acquired bacterial pneumonia and acute bacterial skin/skin structure infections. For newly marketed antimicrobials, few susceptibility (S) test methods are initially available and commercial products (especially automated systems) may not be FDA-approved for more than a year. Interim strategies include testing by agar diffusion methods (disk and Etest) or applying a surrogate test agent that affords high predictive accuracy, especially for MRSA.

Methods: To identify surrogate β -lactams to predict CPT-S, 14,902 USA organisms (3954 SPN; 769 *H. influenzae* [HI]; 8619 *S. aureus* [SA]; 53.6% MRSA); and 1560 indicated Enterobacteriaceae [ENT]) were S tested by CLSI broth microdilution and interpretations, except for CPT (FDA package insert).

Results: For SPN, S to ceftriaxone (CRO) or cefepime (CPM) or amoxicillin/clavulanate (A/C) would accurately predict CPT-S at rates of 99.97, 99.97 and 99.94%, respectively (87.17% of CRO-R pneumococci were CPT-S). HI S to CPT would be predicted at high rates (99.87%) using CRO or CPM or ceftazidime S results. For SA, CPT-S was best predicted by carbapenems (imipenem [IMP] at 99.86%; meropenem [MER] at 99.89%) and IMP- or MER-S or Intermediate (I; MIC, ≤ 8 $\mu\text{g/ml}$) correlated with CPT-S with 99.75 - 99.83% accuracy. CRO was the best surrogate for CPT when testing ENT (95.89%).

Conclusions: CPT can accurately be S tested using a β -lactam surrogate strategy for commercial systems generating MIC results (Vitek[®] 2, Phoenix[™], MicroScan[®]). Accuracy among SPN (99.94-99.97%), HI (99.87%), and SA (99.75-99.89%) was highly acceptable. Among SA 73.99-79.36% of IMP- or MER-R strains remained CPT-S, further minimizing predictive risk and allowing CPT clinical use with local *in vitro* test guidance.

Introduction

Ceftaroline (CPT), the active form of CPT fosamil, is a broad-spectrum cephalosporin with a uniquely high binding affinity for the altered penicillin-binding protein (PBP2a) responsible for methicillin resistance among staphylococci. Unlike other cephalosporins that are inactive against methicillin-resistant *Staphylococcus aureus* (MRSA), ceftaroline has demonstrated *in vitro* potency and clinical success against this important pathogen and has a clinical indication for use in acute bacterial skin and skin structure infections as well as community-acquired bacterial pneumonia (CABP; not MRSA). Therefore, the use of oxacillin and/or ceftazidime test results to predict ceftaroline resistance or susceptibility among other β -lactams does not apply and direct testing of this new cephalosporin would be desirable to predict clinical success per criteria approved in the United States Food and Drug Administration (USA-FDA) product package insert. Those USA-FDA staphylococcal susceptibility criteria are a ceftaroline MIC at ≤ 1 $\mu\text{g/ml}$ and a zone diameter of ≥ 24 mm when using methods published by the Clinical and Laboratory Standards Institute (CLSI). Non-susceptible results (MIC, ≥ 2 $\mu\text{g/ml}$) have not been characterized as either intermediate or resistant due to limited clinical experience with infections caused by staphylococci having those MIC levels.

Similarly, ceftaroline exhibits a high binding to altered PBPs associated with β -lactam MIC elevations in streptococci, particularly *Streptococcus pneumoniae*. This resulting potency advantage compared to ceftriaxone expands the spectrum of ceftaroline against this CABP pathogen and translates to high clinical success rates.

As the vast majority (>80%) of clinical microbiology laboratories do not use reference/standardized CLSI methods, alternative strategies for testing newly released antimicrobials must be developed due to long term delays in the development and USA-FDA approval of commercial susceptibility testing products (Vitek[®], Vitek[®] 2, BD Phoenix[™], MicroScan[®], Sensititre[®]). As some antimicrobials may present immediate therapeutic advantages, one strategy is to test a surrogate agent (usually in the same class) as a predictor of susceptibility and/or resistance. This testing option has been most recently applied to doripenem, but has also been used for other β -lactams (cefotetan, cefpodoxime), and by the CLSI in Table 1 of document M100-S22. The most difficult obstacle for ceftaroline has been to select an appropriate antimicrobial class (β -lactam) agent when testing MRSA strains, where no other commercially available β -lactam has demonstrated clear *in vitro* and clinical utility. However, some carbapenems have shown measurable potencies versus MRSA that may be usable, as would advanced-spectrum cephalosporins (ceftriaxone, ceftazidime or cefepime) when testing *Streptococcus* spp. This study investigated the optimal use of candidate β -lactams as surrogate predictors of ceftaroline activity/susceptibility, allowing the earliest guided clinical use in medical centers having FDA-approved commercial susceptibility testing systems reporting quantitative MIC values or category interpretations using the CLSI and USA-FDA breakpoint criteria.

Methods

The organisms tested in the Assessing Worldwide Antimicrobial Resistance Evaluation (AWARE) Surveillance Program from 2008-2010 were analyzed to select a surrogate marker agent for ceftaroline. The species selected were: 3954 *S. pneumoniae* (2008-2010); 769 *Haemophilus influenzae* (2010 only); 8619 *S. aureus* (2008-2010; 53.6% MRSA); and 1560 indicated Enterobacteriaceae (2010 only, *E. coli* and *Klebsiella* spp), for a total of 14,902 strains, all tested by the CLSI M07-A9 method in a GLP facility (JMI Laboratories, North Liberty, Iowa). Concurrent quality control (QC) strains *E. coli* ATCC 25922, *S. aureus* ATCC 25923 and 29213, *H. influenzae* ATCC 49247 and *S. pneumoniae* ATCC 49619 were used, and all QC results were within published ranges for ceftaroline and candidate surrogate β -lactams.

Analysis focused on the identification of surrogate agents to predict ceftaroline susceptibility minimizing, where possible, false-susceptibility to $\leq 1.5\%$ and false-intermediate rates to $\leq 5\%$. Comparisons used published breakpoint criteria for each agent. With the exception of testing the Enterobacteriaceae, ceftaroline only has susceptible and non-susceptible criteria precluding total cross-resistance calculations. By using the most potent β -lactam surrogate agents against each pathogen, the interpretive error was found to be far below the listed target accuracy limits above (see Tables and Figures).

Results

Surrogate Candidates for Staphylococci

– Carbapenems (imipenem and meropenem) demonstrated measurable activities against the tested MRSA (53.6% of *S. aureus* sample), having MIC values ranging from ≤ 0.12 to >8 $\mu\text{g/ml}$ for both agents. This high carbapenem activity against MRSA and all *S. aureus* made them unique surrogate candidates among β -lactams and demonstrated a linear relationship to ceftaroline MIC results (Figures 1 and 2)

– For imipenem used as a surrogate for ceftaroline susceptibility (Figure 1), an MIC of ≤ 4 $\mu\text{g/ml}$ to imipenem was 99.86% accurate in predicting a susceptible ceftaroline MIC value (≤ 1 $\mu\text{g/ml}$). By also adding imipenem MIC results at 8 $\mu\text{g/ml}$ (eg. ≤ 8 $\mu\text{g/ml}$ for imipenem = ceftaroline-susceptible), the accuracy rate only decreased to 99.75%. The accuracy of ceftaroline-susceptible values predicted by meropenem MIC results at ≤ 4 and ≤ 8 $\mu\text{g/ml}$ was 99.89 and 99.83%, respectively (Figure 2)

– Among imipenem- and meropenem-non-susceptible *S. aureus* tested (MICs, ≥ 8 $\mu\text{g/ml}$), 73.99 and 79.36% of isolates (all MRSA) remained susceptible to ceftaroline

Surrogate Candidates for *S. pneumoniae*

– Only the most active commercially available β -lactams would be possible surrogates for ceftaroline susceptibility and their accuracy is shown in Table 1. The accuracy rates of ceftriaxone (99.97%) and cefepime (99.97%) were best to predict ceftaroline-susceptible *S. pneumoniae*; however, amoxicillin/clavulanate (99.94%) could also be used with confidence. Note that utilizing amoxicillin/clavulanate susceptibility results significantly underestimated coverage of ceftaroline when compared to using ceftriaxone (Figure 3 and Table 1)

– Also no ceftriaxone-intermediate strains (336 isolates) were ceftaroline-non-susceptible (Table 1); therefore a ceftriaxone MIC at ≤ 2 $\mu\text{g/ml}$ predicted ceftaroline susceptibility with an accuracy of nearly 100.0%

Surrogate Candidates for *H. influenzae*

– All three agents (ceftriaxone, cefepime, and ceftazidime) performed well (99.87% accuracy) and can be used with confidence as surrogate agents for ceftaroline susceptibility versus *H. influenzae*; even with the very low USA-FDA ceftaroline breakpoint at ≤ 0.12 $\mu\text{g/ml}$, compared to the other cephalosporins (≤ 2 $\mu\text{g/ml}$ as susceptible; see Table 2)

Surrogate Candidates for Indicated Enterobacteriaceae

– A total of 1560 *E. coli* and *Klebsiella* spp. isolated in 2010 were used to assess the accuracy of using ceftriaxone (Figure 4) and other cephalosporins to predict ceftaroline susceptibility at ≤ 0.5 $\mu\text{g/ml}$. Using ceftriaxone to predict ceftaroline-susceptible strains showed a 95.89% accuracy with 2.09 and 2.02% minor and very major errors, respectively. If the ceftaroline-susceptible breakpoint was adjusted to ≤ 1 $\mu\text{g/ml}$ (2 $\mu\text{g/ml}$ as intermediate), like that used for *S. aureus*, the accuracy rate would be 97.98% for ceftriaxone with only 1.37% very major error (acceptable), see Figure 4.

Table 1. Susceptibility Category Comparisons of Ceftaroline with Three Broad-spectrum β -lactams (Ceftriaxone, Cefepime, Amoxicillin/clavulanate) when Tested against 3954 *S. pneumoniae* from the USA

Surrogate candidate	CLSI category (MIC)	No. ceftaroline MICs by category ($\mu\text{g/ml}$):		
		Susceptible (≤ 0.25)	(0.5)	(≥ 1)
Ceftriaxone	Susceptible (≤ 1)	3532 ^a	1 ^a	0
	Intermediate (2)	336	0	0
	Resistant (≥ 4)	31	54	0
Cefepime	Susceptible (≤ 1)	3537 ^b	1 ^b	0
	Intermediate (2)	354	29	0
	Resistant (≥ 4)	8	25	0
Amoxicillin/clavulanate	Susceptible (≤ 1)	3283 ^c	2 ^c	0
	Intermediate (2)	107	4	0
	Resistant (≥ 4)	509	49	0

a. Accuracy of ceftriaxone susceptibility results to predict ceftaroline susceptibility at 3532/3533 (99.97%); among the 421 ceftriaxone-non-susceptible pneumococcal strains, 87.17% were ceftaroline-susceptible.
b. Accuracy at 99.97% (3537/3538).
c. Accuracy at 99.94% (3283/3285).

Table 2. Susceptibility Category Comparisons of Ceftaroline with Three Cephalosporins (Ceftriaxone, Cefepime, Ceftazidime) when Tested against 769 *H. influenzae* from the USA

Surrogate candidate	CLSI category (MIC)	No. ceftaroline MICs by category ($\mu\text{g/ml}$):	
		Susceptible (≤ 0.12)	Non-susceptible (≥ 0.25)
Ceftriaxone	Susceptible (≤ 2)	768 ^a	1 ^a
	Non-susceptible (≥ 4)	0	0
Cefepime	Susceptible (≤ 2)	768 ^a	1 ^a
	Non-susceptible (≥ 4)	0	0
Ceftazidime	Susceptible (≤ 2)	768 ^a	1 ^a
	Non-susceptible (≥ 4)	0	0

a. Accuracy of ceftriaxone, cefepime and ceftazidime susceptibility results to predict ceftaroline susceptibility was at 768/769 (99.87%).

Figure 1. Scattergram of Ceftaroline MIC Values Compared to Imipenem MIC Results when Tested against 8619 *S. aureus* (4624 or 53.6% were MRSA) from the USA. Solid Bolded Horizontal (USA-FDA) and Vertical (CLSI) Lines Indicate Breakpoints for Each Agent




Figure 2. Scattergram of Ceftaroline MIC Values Compared to Meropenem MIC Results when Tested against 8619 *S. aureus* (4624 or 53.6% were MRSA) from the USA. Solid Bolded Horizontal (USA-FDA) and Vertical (CLSI) Lines Indicate Breakpoints for Each Agent




Figure 3. Scattergram of Ceftaroline MIC Values Compared to Ceftriaxone MICs when Testing 3954 *S. pneumoniae* Isolates from the USA. Only 1.39% of Pneumococci were Ceftaroline Non-susceptible and No MIC was Observed at >0.5 $\mu\text{g/ml}$. Solid Bolded Horizontal (USA-FDA) and Vertical (CLSI) Lines Indicate Breakpoints for Each Agent





Figure 4. Scattergram Comparing Ceftaroline and Ceftriaxone MIC Values when Testing 1560 Indicated Enterobacteriaceae (*E. coli* and *Klebsiella* spp.). Two Ceftaroline Breakpoints are Shown, 1.) USA-FDA Product Package Insert (Solid Bolded Horizontal Lines) and 2.) ≤ 1 $\mu\text{g/ml}$ as Susceptible, Like that Applied to Staphylococci (Broken Horizontal Line)



Conclusions

• This interim ceftaroline susceptibility testing strategy of using potent β -lactams as surrogate agents was examined using 14,902 recent clinical isolates tested by reference CLSI broth microdilution methods and published breakpoint concentrations (CLSI and USA-FDA)

• High-level accuracy for predicting ceftaroline susceptibility was observed as follows: for *S. aureus* use either imipenem or meropenem at a susceptible MIC of ≤ 8 $\mu\text{g/ml}$ (99.75-99.83% accurate); for *S. pneumoniae* use ceftriaxone at a susceptible MIC of ≤ 2 $\mu\text{g/ml}$ or at ≤ 1 $\mu\text{g/ml}$ ($\geq 99.97\%$ accurate); for *H. influenzae* use ceftriaxone or cefepime or ceftazidime at a susceptible MIC of ≤ 2 $\mu\text{g/ml}$ (99.87% accurate); and for indicated *E. coli* and *Klebsiella* spp. use ceftriaxone at a susceptible MIC of ≤ 1 $\mu\text{g/ml}$ (95.89% accurate)

• Like previous uses of surrogate agents within the β -lactam antimicrobial class, these qualified uses of carbapenems (*S. aureus*) and "third-or fourth-generation" cephalosporins (*S. pneumoniae*, *H. influenzae*, and indicated Enterobacteriaceae) for ceftaroline testing offers immediate applications to existing susceptibility testing results generated by commercial systems producing quantitative (MICs) or categorical/qualitative results.

References

- Barry AL, Jones RN (1987). Cross susceptibility and absence of cross resistance to cefotetan and ceftioxin. *J Clin Microbiol* 25: 1570-1571.
- Clinical and Laboratory Standards Institute (2012). *M07-A9. Methods for dilution antimicrobial susceptibility tests for bacteria that grow aerobically; approved standard: ninth edition*. Wayne, PA: CLSI.
- Clinical and Laboratory Standards Institute (2012). *M100-S22. Performance standards for antimicrobial susceptibility testing: 22nd informational supplement*. Wayne, PA: CLSI.
- Corey GR, Wilcox M, Talbot GH, Friedland HD, Baculik T, Withereil GW, Critchley I, Das AF, Thye D (2010). Integrated analysis of CANVAS 1 and 2: Phase 3, multicenter, randomized, double-blind studies to evaluate the safety and efficacy of ceftaroline versus vancomycin plus aztreonam in complicated skin and skin-structure infection. *Clin Infect Dis* 51: 641-650.
- File TM, Jr., Low DE, Eckburg PB, Talbot GH, Friedland HD, Lee J, Llorens L, Critchley I, Thye D (2010). Integrated analysis of FOCUS 1 and FOCUS 2: Randomized, double-blind, multicenter phase 3 trials of the efficacy and safety of ceftaroline fosamil versus ceftriaxone in patients with community-acquired pneumonia. *Clin Infect Dis* 51: 1395-1405.
- Jones RN, Farrell DJ, Mendes RE, Sader HS (2011). Comparative ceftaroline activity tested against pathogens associated with community-acquired pneumonia: Results from an international surveillance study. *J Antimicrob Chemother* 66 Suppl 3: iii69-iii80.
- Jones RN, Mendes RE, Sader HS (2010). Ceftaroline activity against pathogens associated with complicated skin and skin structure infections: Results from an international surveillance study. *J Antimicrob Chemother* 65 Suppl 4: iv17-iv31.
- Jones RN, Sader HS, Fritsche TR, Janecek MJ (2007). Selection of a surrogate beta-lactam testing agent for initial susceptibility testing of doripenem, a new carbapenem. *Diagn Microbiol Infect Dis* 59: 467-72.
- Jones RN, Zurenko GE (1993). Prediction of bacterial susceptibility to cefpodoxime by using the ceftazidime minimum inhibitory concentration result. *Diagn Microbiol Infect Dis* 17: 313-316.
- Teflaro Package Insert (2010). Available at http://www.frx.com/pi/Teflaro_pi.pdf. Accessed April 14, 2011.
- Zhanell GG, Snieszek G, Schweizer F, Zelenitsky S, Lagace-Wiens PR, Rubinstein E, Gin AS, Hoban DJ, Karlowsky JA (2009). Ceftaroline: A novel broad-spectrum cephalosporin with activity against methicillin-resistant *Staphylococcus aureus*. *Drugs* 69: 809-831.

Acknowledgments

RNJ, RKF, HSS and MGS are employees of JMI Laboratories that coordinates the AWARE ceftaroline surveillance program for Cerexa/Forest Laboratories Inc. JMI Laboratories receives funding from various pharmaceutical industry sources for *in vitro* evaluations/surveillance of β -lactam agents including those analyzed in this publication. RNJ is the guarantor of the data. The co-authors also acknowledge the following JMI employees for support during the analysis and preparation of this manuscript: S. Benning, P. Clark, D. Farrell, A. Fuhrmeister, and G. Moet.

This study was supported by Forest Laboratories, Inc. Forest Laboratories, Inc. was involved in the study design and decision to present these results. Forest Laboratories, Inc. had no involvement in the analysis, collection, and interpretation of data.

Scientific Therapeutics Information, Inc., provided editorial assistance which was funded by Forest Research Institute, Inc.