

Pilot Studies to Guide CLSI Re-evaluations of *S. pneumoniae* and β -haemolytic Streptococci Disk Diffusion Criteria for Commonly Used Penicillins and Cephalosporins

RN JONES¹, PR RHOMBERG¹, JD TURNIDGE²

¹JMI Laboratories, North Liberty, IA, USA; ²Woman's & Children's Hospital, Adelaide, Australia

ABSTRACT

Background: CLSI β -lactam (β -L) disk diffusion (DD; Kirby-Bauer) methods have no correlate zone for the old or new susceptible (S) breakpoints for SPN, and β -haemolytic streptococci (BHS)-S criteria have been guided by penicillin (PEN)-S results. With the recently modified PEN-S criteria (≤ 2 μ g/ml) for high-dose pneumococcal (SPN) pneumonia therapy, the DD test was re-examined using alternative DD contents for 3 agents and X-S for 12 β -L drugs.

Methods: 153 SPN and 50 BHS were selected from year 2008 isolates to represent all relevant wildtype and S or R subsets. PEN (1-, 2-, 10-U), cefotaxime (CTX; 5-, 30- μ g), ceftriaxone (CRO; 5-, 30- μ g) were tested with several candidate disk contents, and amoxicillin/clavulanate (AUG) and cefepime (FEP) were also tested. BHS were processed by X-S analyses to PEN by MIC and zone criteria (12 β -lactams). All tests used CLSI M02-A10, M07-A8 (2009) methods and M100-S19 criteria.

Results: BHS X-S results indicate PEN (MIC, ≤ 0.12 μ g) continues to predict high activity (MIC, ≤ 0.25 μ g/ml) for AUG, CTX, CRO, FEP, cephalothin and cefuroxime, but *not* cephalaxin (MIC range to 4 μ g/ml), cefadroxil (2 μ g/ml), cefaclor (1 μ g/ml) and ceftazidime (2 μ g/ml). SPN PEN X-S analyses at ≤ 2 μ g/ml favorably correlated with CTX, CRO, FEP and AUG activity. DD testing could be applied with acceptably low error rates for 5 commonly used β -lactams (Table), and X-S was without serious false-S or false-resistance error (86.9-92.2% absolute categorical agreement).

Table: Examples of documented acceptable DD-S breakpoint criteria when testing SPN against five popular antimicrobials.

Agent	Disk conc. (U or μ g)	S Criteria		Error rates (%) ^a			Acceptable? (%)
		MIC (μ g/ml)	DD zone (mm)	VM	Ma	Mi	
PEN	1	≤ 0.06	≥ 24	1.3	0.0	-	Yes (98.7)
	2	≤ 0.06	≥ 26	0.6	2.6	-	Yes (96.8)
	10	≤ 0.06	≥ 32	0.6	2.0	-	Yes (97.4)
	10	≤ 2	≥ 19	0.0	0.0	7.2	Yes (92.8)
CRO	5	≤ 1	≥ 19	0.0	0.0	3.3	Yes (96.7)
	30	≤ 1	≥ 26	0.0	0.0	4.6	Yes (95.4)
	30	≤ 1	≥ 19	0.0	0.0	4.6	Yes (95.4)
CTX	5	≤ 1	≥ 27	0.0	0.0	5.2	Yes (94.8)
	30	≤ 1	≥ 26	0.0	0.0	7.2	Yes (92.8)
AUG	30	≤ 2	≥ 22	0.0	0.0	5.2	Yes (94.8)

a. VM = very major (false-S), Ma = major (false-R) and Mi = minor.

Conclusions: DD testing for contemporary SPN & BHS appear practical and accurate for most widely used β -L's, demonstrating acceptable performance and potential X-S criteria use (PEN surrogate). CLSI should re-establish simple, cost effective DD criteria for these prevalent pathogens.

INTRODUCTION

With the adoption of the Clinical and Laboratory Standards [CLSI; formerly the National Committee for Clinical Laboratory Standards (NCCLS)] breakpoint tables in the 1970's, the recommendation for testing one antimicrobial agent to represent the susceptibility for other similar agents has been widely accepted (class-disk concept). One such instance is the testing of a penicillin disk (10-U) against *S. pneumoniae* and β -haemolytic streptococcal isolates to predict the susceptibility to other β -lactam agents including β -lactam/ β -lactamase inhibitor combinations and carbapenems. Similarly the cephalothin (30- μ g) or cefazolin

(30- μ g) disks have been used to predict susceptibility (not resistance) to other cephalosporin compounds, especially orally delivered agents.

With the recently adopted breakpoints for penicillin susceptibility of ≤ 2 μ g/ml (high-dose regimens; 12 million units/day) and resistance at ≥ 8 μ g/ml against *S. pneumoniae*, the disk diffusion test should be re-evaluated to confirm the cross-susceptibility with other β -lactam agents since no disk zone diameter breakpoints are found in the current CLSI M100-S19 tables.

MATERIALS AND METHODS

Bacterial Isolates. A collection of 50 β -haemolytic streptococci (primarily *S. pyogenes*, 80%) and 153 *S. pneumoniae* were selected from the 2008 SENTRY Antimicrobial Resistance Surveillance study based on their susceptibility to penicillin. *S. pneumoniae* isolates were evenly distributed according to non-meningitis, oral penicillin breakpoints; susceptible at ≤ 0.06 μ g/ml, intermediate at 0.12-1 μ g/ml and resistant at ≥ 2 μ g/ml.

Susceptibility Testing. Frozen reference broth microdilution panels with Mueller-Hinton broth supplemented with 2-5% lysed horse blood were produced by TREK Diagnostics (Cleveland, Ohio, USA) containing penicillin (PEN), ampicillin, amoxicillin/clavulanate, cefaclor, cephalaxin, cefadroxil, cephalothin, cefuroxime, ceftazidime, ceftriaxone, cefotaxime, and cefepime were inoculated according to CLSI standard methods (M07-A8, 2009). Concurrently, the CLSI disk diffusion (DD) method (M02-A10, 2009) with 150mm diameter Mueller-Hinton agar plates supplemented with 5% sheep blood were used for testing a total of nine disks including penicillin (1-, 2-, 10-U), amoxicillin/clavulanate (30- μ g), ceftriaxone (5-, 30- μ g), cefotaxime (5-, 30- μ g) and cefepime (30- μ g).

RESULTS

- Against BHS, the PEN MIC test (susceptible breakpoint at ≤ 0.12 μ g/ml) accurately predicted high activity (MIC, ≤ 0.25 μ g/ml) for six of the β -lactam agents tested, but *not* for cephalaxin (MIC values ranging to 4 μ g/ml), cefadroxil (range to 2 μ g/ml), ceftazidime (range to 2 μ g/ml) and cefaclor (range to 1 μ g/ml; Table 1).
- The MIC/DD zone diameter correlations for the seven β -lactam agents tested against BHS all showed MIC values ≤ 0.25 μ g/ml and zone diameter results ≥ 21 mm (Data not shown).

- The new PEN-susceptible MIC breakpoint of ≤ 2 μ g/ml for high-dose *S. pneumoniae* therapy showed excellent predictive cross-susceptibility with cefotaxime (minor error rate only, 7.8%), cefepime (8.5%) and ceftriaxone (9.2%). Elevated unacceptable error rates were observed for ceftazidime (12.4%), amoxicillin/clavulanate (13.1%), ampicillin (13.8%), cefuroxime axetil (40.0%) and cefaclor (50.9%; Table 2).

- For the *S. pneumoniae* isolates, the MIC/DD correlation for the five β -lactam agents tested showed an acceptable susceptible predictive ability with the tested lower disk concentrations for PEN, ceftriaxone and cefotaxime, each generally performing better than current CLSI-recommended disk drug contents. Accuracy rates ranged from 92.8 to 98.7% (Table 3).

- Excellent intermethod correlations of MIC/DD zone diameter results using regression analysis for the *S. pneumoniae* tests were observed. Solid lines at proposed susceptible and resistant MIC and DD breakpoints show the limited number of strains producing error rates, only 1.3-7.2% (Table 3; Figures 1 and 2).

Table 1. Penicillin cross-resistance data for 50 β -haemolytic streptococci tested using CLSI broth microdilution methods against 11 β -lactam agents.

Compound Agent	MIC (μ g/ml)	Penicillin MIC (μ g/ml):							% susceptible ^b
		≤ 0.008	0.015	0.03	0.06	0.12	0.25		
Ampicillin	≤ 0.5	41	3	5	1	-	-	-	100.0 ^c
Amoxicillin/clavulanate	≤ 0.12	41	3	5	1	-	-	-	100.0 ^c
Cephalothin	≤ 0.12	41	2	5	-	-	-	-	100.0 ^c
Cephalaxin	≤ 0.25	33	-	-	-	-	-	-	^c
Cefadroxil	≤ 0.25	41	3	-	-	-	-	-	^c
Cefaclor	≤ 0.12	41	2	-	-	-	-	-	100.0 ^b
Cefuroxime	≤ 0.12	41	3	5	1	-	-	-	100.0 ^b
Cefotaxime	≤ 0.06	41	3	4	1	-	-	-	100.0
Ceftriaxone	≤ 0.06	41	3	4	1	-	-	-	100.0
Ceftazidime	≤ 0.12	41	3	3	1	-	-	-	^c
Cefepime	≤ 0.06	41	3	4	1	-	-	-	100.0

a. CLSI M100-S19 breakpoint criteria were used.
b. Susceptibility rate based on penicillin MIC per CLSI recommendations.
c. No interpretive criteria provided by CLSI.
d. Boxed values may indicate non-susceptible MIC values.

Table 2. Penicillin cross-susceptibility data for 153 *S. pneumoniae* tested using CLSI methods against 11 β -lactam agents.

Compound Agent	MIC (μ g/ml)	Penicillin MIC (μ g/ml):											% susc. ^b						
		≤ 0.008	0.015	0.03	0.06	0.12	0.25	0.5	1	2	4 ^a	8		16	>16				
Ampicillin	≤ 0.5	16	25	6	6	14	13	10	1	-	-	-	-	-	-	-	-	-	34.6 ^c
Amox/clav	≤ 0.12	16	25	6	6	14	9	4	-	-	-	-	-	-	-	-	-	-	86.7
Cefaclor	≤ 0.25	4	3	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	69.9
Cefuroxime	≤ 0.2	15	25	4	6	5	1	1	1	1	-	-	-	-	-	-	-	-	60.1
Cefotaxime	≤ 0.06	15	25	5	6	10	3	2	-	-	-	-	-	-	-	-	-	-	69.9
Ceftriaxone	≤ 0.06	16	25	6	6	12	4	1	-	-	-	-	-	-	-	-	-	-	60.1
Cefepime	≤ 0.06	16	25	6	6	9	2	-	-	-	-	-	-	-	-	-	-	-	69.9
Cefadroxil	≤ 0.25	16	16	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	^c
Ceftazidime	≤ 0.12	16	16	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	^c

a. CLSI M100-S19 breakpoint criteria were used.
b. Solid vertical and horizontal lines represent susceptible and resistant MIC breakpoint criteria.
c. Based on oral penicillin breakpoint at ≤ 0.06 μ g/ml.
d. No interpretive criteria provided by CLSI.

Table 3. Disk diffusion susceptibility breakpoint criteria compared to MIC breakpoints when testing *S. pneumoniae* isolates against five β -lactam antimicrobials.

Compound Agent	Disk conc. (U or μ g)	Susceptibility Criteria		Intermethod Error rates (%) ^a			Acceptable? (% Agreement)	Associated Figure
		MIC (μ g/ml)	DD zone (mm)	VM	Ma	Mi		
Penicillin	1	≤ 0.06	≥ 24	1.3	0.0	-	Yes (98.7)	1
	2	≤ 0.06	≥ 26	0.6	2.6	-	Yes (96.8)	-
	10	≤ 0.06	≥ 32	0.6	2.0	-	Yes (97.4)	-
Ceftriaxone	5	≤ 1	≥ 19	0.0	0.0	7.2	Yes (92.8)	1
	30	≤ 1	≥ 26	0.0	0.0	3.3	Yes (96.7)	2
	30	≤ 1	≥ 19	0.0	0.0	4.6	Yes (95.4)	2
Cefotaxime	5	≤ 1	≥ 26	0.0	0.0	4.6	Yes (95.4)	-
	30	≤ 1	≥ 27	0.0	0.0	5.2	Yes (94.8)	-
	30	≤ 1	≥ 26	0.0	0.0	7.2	Yes (92.8)	2
Amoxicillin/clavulanate	30	≤ 2	≥ 23	0.0	0.0	3.9	Yes (96.1)	2

a. VM = very major (false-S), Ma = major (false-R) and Mi = minor.

Figure 1. Penicillin MIC versus disk diffusion zone diameter (1- and 10-U disks) for 153 *S. pneumoniae* tested using CLSI methods.





Figure 2. MIC versus disk diffusion zone diameter results for other β -lactam agents against *S. pneumoniae* isolates tested using CLSI methods.



CONCLUSIONS

- The CLSI disk diffusion method to test β -lactam agents against contemporary clinical *S. pneumoniae* and β -haemolytic streptococci isolates can be used against contemporary strains, and the penicillin disk demonstrates acceptable cross-susceptibility predictive accuracy for other selected β -lactam agents.

- The use of lower drug disk concentrations for penicillin (1-U), ceftriaxone (5- μ g) and cefotaxime (5- μ g) could be used to increase the overall agreement between the broth microdilution and disk diffusion methods. Breakpoints for these DD tests are proposed.

- CLSI should consider re-establishing disk diffusion susceptibility criteria for testing *S. pneumoniae* and β -haemolytic streptococci against some β -lactam agents (penicillins \pm enzyme inhibitors, cefepime, cefotaxime and ceftriaxone), as a simple cost effective and accurate method.

ACKNOWLEDGEMENTS

The authors would like to thank Becton/Dickinson, Sparks, MD, USA; Mast Group Ltd., Merseyside, UK, and TREK Diagnostics, Cleveland, OH, USA for providing the agar test media, antimicrobial disks and broth microdilution reference GMP panels for this study. Technical and editorial assistance was provided by the following individuals; M.G. Stilwell, M. J. Janeczek, D.J. Biedenbach and A. Small.

REFERENCES

- Clinical and Laboratory Standards Institute (2009). M02-A10. Performance standards for antimicrobial disk susceptibility tests; approved standard - tenth edition. Wayne, PA:CLSI.
- Clinical and Laboratory Standards Institute (2009). M07-A8. Methods for dilution antimicrobial susceptibility tests for bacteria that grow aerobically; approved standard - eighth edition. Wayne, PA:CLSI.
- Clinical and Laboratory Standards Institute (2009). M100-S19. Performance standards for antimicrobial susceptibility testing. 19th informational supplement. Wayne, PA:CLSI.
- Weinstein MP, Klugman KP, Jones RN (2009). Rationale for revised penicillin susceptibility breakpoints versus *Streptococcus pneumoniae*: coping with antimicrobial susceptibility in an era of resistance. *Clin Infect Dis* 48: 1596-1600.